



CryoTrans series

Introduction

CryoTrans series is designed for storage and short-distance transportation of small amount of liquid nitrogen. It is equipped with rubber protection rings and prefixed bottom pad for safety. Stainless steel roller base is optional for convenient transportation. CryoTrans series is widely used in animal husbandry and laboratories.

Key Features

- 1** Strong, lightweight aluminum construction
- 2** Low liquid nitrogen evaporation
- 3** Special liquid nitrogen transportation design
- 4** CE Certificate
- 5** 5-year vacuum warranty
- 6** Equipped with rubber protection rings & prefixed bottom pad

Important Accessories

- 1. Liquid nitrogen level ruler
- 2. Liquid nitrogen dispenser
- 3. Roller base.



Technical Parameters

Model	CryoTrans 3	CryoTrans 6	CryoTrans 10	CryoTrans 20	CryoTrans 30	CryoTrans 35	CryoTrans 50
Performance							
Capacity (L)	3	6	10	20	30	35	50
Neck Diameter (mm)	50	50	50	50	50	50	50
Static Evaporation Rate (L/day)	0.12	0.12	0.12	0.12	0.12	0.12	0.23
Unit Dimensions							
Overall Height (mm)	435	482	552	672	706	750	811
External Diameter (mm)	223	300	300	394	462	462	462
Weight Empty (KG)	3.1	4.8	5.9	9.5	12.9	14.2	15.4
Weight Full (KG)	5.56	9.72	14.1	25.9	37.5	42.9	56.4