JENCO









Advance Portable Meter

Dissolved Oxygen, Temp.

The Jenco 9250 and 9251 microprocessor based dissolved oxygen meters give you laboratory grade measurements in an easy to use, portable package:

Easy and Accurate Measurements

Clear LCD displays D.O. as percent air saturation or D.O. as parts per thousand, simultaneously with temperature.

Manual or automatic temperature compensated (ATC) D.O. readings.

User selectable salinity and pressure compensated D.O. readings.

Works with accurate and reliable polarographic D.O. probes.

Convenient one point calibration for D.O. Calibration data is stored in memory and is ready for use on power up.

Automatic lock mode and end point sensing are available to lock in stable readings.

Intuitive operation with membrane keypad with tactile and audio feedback.

Inclined back provides ideal viewing angle for benchtop use.

Carefree Maintenance

Watertight case is splash resistant and made with durable plastic intended for rugged use. Protective covers for connectors are included.

Standard Inputs and Outputs

RS-232 computer interface with software (3250 only).

8 pin DIN connector for D.O. & thermistor.

AC adaptor or battery power.

9250 / 9251 Product Specifications

Specifications

D.O. (ppm/air saturated)

Range

0.00 to 20.00 ppm / 0 to 199.9 %

Resolution

0.01 ppm / 0.1 %

± 0.2 % full scale

Temp -6.0 to 46.0 °C 0.1 °C $\pm 0.1 \text{ °C} \pm 1 \text{ LSD}$

General

D.O.:

Salinity Compensation Manual 0.0 to 40.0 ppt
Pressure Compensation Manual 640 to 1100 mBar
Temperature Compensation Auto -10.0 to 50.0 °C
Sensor Type Polarographic

Temperature:

Sensor Type Thermistor 2.252 k ohm

Other:

Power Supply 6 AAA Batteries or 115/230 V AC adapter

Battery Life ~ 25 hours

Calibration Backup Yes

Audio Feedback All Touch Keys

End Point Sensing & Hold Yes

LCD Character Size (Main, Temp) 15 mm, 8.7 mm high LCD

Ambient Temp. Range 0 to 50 °C RS-232 communication Bate Rate 9600

Dimensions (W x D x H) 86 mm x 226 mm x 72 mm
Weight 410 grams (Including batteries)

Accessories

900YDOL-10 Watertight D.O. probe w/ 10ft cable

900YC-25 25ft cable

9005YDOK-L D.O. powder, sandpaper and 6 caps

6007D115 115 V AC adapter 6007D230 230 V AC adapter

